


# TECHNICAL NOTE

TN0007: Optimiser™ Microplate System (ELISA) Setup Guide  
on the PerkinElmer VICTOR™ X Multilabel Microplate Reader



*Better Immunoassays through  
Innovative Microfluidics*

## READER SETUP WITH PERKINELMER VICTOR™ X MULTILABEL PLATE READER :

The PerkinElmer VICTOR™ X Multilabel Plate Reader with fluorescence function has been tested for compatibility with Siloam's Optimiser™ microplate System in ELISA assay with OptiGlow™ chemifluorescent substrate. Please refer to the Optimiser™ Technology page on Siloam's website for more details on the principles behind the Optimiser™ microplate platform. For more detailed information of Infinite® readers or software, please contact PerkinElmer technical support.

### Recommended Optics

	Product number	Wavelength
<b>Excitation</b>	1420-5680	531/25 nm or similar
<b>Emission</b>	1420-5860	590/10 nm or similar


*Note: Please refer to Spectra of OptiGlow™ substrate in General Reader Setup Guide under Optimiser™ Technology page for selecting filters.*

### Instrument Setup

Turn on the plate reader, and open up PerkinElmer 2030 Manager. Please refer to reader instruction manual for detail operation.


#### Step 1- Edit Fluorometry Label Properties

In panel of "Labels", choose "**Fluorometry**", and add a new label file.



In the new Label, select or input the properties as listed below:


- CW-Lamp Energy:** 1000
- CW-Lamp Control:** Stabilized Energy
- CW-Lamp Filter:** Select the appropriate excitation filter
- Excitation Aperture:** Normal
- Emission Filter:** Select the appropriate emission filter
- Emission Aperture:** Normal
- Counter position:** Top
- Counting Time:** 0.1



*Note: For VICTOR X reader, the absolute value of fluorescence intensity is relative to the setting of CW-Lamp Energy and Counting Times. But it won't affect the signal/noise ratio and the assay performance.*

## Step 2- Edit the Protocol

Open the “Protocol Editor”, select “Generic 8x12 size plate” for Plate type. Select “Standard” for Measurement height.



### Step 3- Verification with Positive Control Solution

#### **Material:**

- 1) Siloam's SAV-HRP (cat# OMR-HRP). *Note: Any HRP conjugate with concentration greater than 1 µg/mL can be used for this test with following experimental protocol.*
- 2) Siloam's OptiGlow™ substrate.
- 3) One Optimiser™ Microplate.

#### **Positive Control Preparation and Loading:**

In a clean plastic tube, add 50 µL of OptiGlow™-A, 50 µL of OptiGlow™-B, 1 µL of OptiGlow™-C, and 1 µL of supplied SAV-HRP stock solution, mix well, and wait for 2 minutes. The substrate will be fully developed and stable for hours. Load 4 µL of mixture into one well (*well A1 is used in this instruction*) of Optimiser™ microplate and wait until the well is empty (do not use pad/holder)

#### **Signal Reading:**

Follow the setting in step 1 and 2, read well A1 and A2. Below is a typical data\* from a VICTOR X reader.

	A1	A2
RFU	51283	190

*\*The actual value of the reading may vary between readers.*

**Technical Assistance:** If you require assistance, please contact Siloam Biosciences, Inc. Technical Support at +1 (513) 429-2976 or [techsupport@siloambio.com](mailto:techsupport@siloambio.com).

Additional technical assistance is available under the Technical Support tab on the Siloam Biosciences web site (<http://siloambio.com/>).

- Using Optimiser™ Immunoassay Microplate Video
- Optimiser™ User's Guide
- Reader Settings
- Quick Reference Guide
- Frequently Asked Questions
- Application Notes

Two additional videos appear under the Technology tab of the web site.

- Optimiser™ Principles of Operation
- Running an Assay with Optimiser™



***Better Immunoassays Through Innovative Microfluidics***

***Siloam Biosciences, Inc.***

413 Northland Blvd., Cincinnati, OH 45240  
USA

Phone: +1 (513) 429-2976

Fax: +1 (513) 429-2946

<http://www.siloambio.com>